
SEPTEMBER 2023 MONTHLY HBOT WEBINAR

Clean Chambers are Safe Chambers

PRESENTED BY THE ACMH PROGRAM





Maintain and Clean Chambers

- Maintaining clean chambers protects their longevity and the health of your patients.
- The chambers should be cleaned with approved cleaners weekly or between patients if the patient has contact precautions due to infection, such as MRSA.
- Always clean the chambers sooner if they become visibly soiled with anything including fingerprints on the acrylic.





How to Clean Chambers

- Use a cotton cloth when wiping the acrylic, such as a pillowcase, and do not clean the acrylic in a circular motion.
- When cleaning the acrylic be sure to check for crazing or scratches, if you notice any acrylic damage, notify your hyperbaric service technician and national safety director right away.





How to Clean Chambers

- Mattresses must be cleaned between each patient with non-alcohol wipes and fresh, clean linens.
- The muffler should be removed from the chamber weekly and cleaned with soapy water and a wire brush, if necessary.





Area around the Chambers

- Make sure that the area surrounding the chamber is always kept clean and free of dust.
- If you notice an accumulation of dust in the door lock of the chamber, attempt to clean the dust away, but let the hyperbaric service technician know if the integrity of the door lock becomes compromised.
- Use a Swiffer or another duster to clean the outside of the chamber including between the bottom of the acrylic and the base of the chamber.



Approved Cleaning Products

- Approved cleaning products can be found at
 - perrybaromedical.com
 - sechristusa.com





Hand Sanitizers

 Hand sanitizers can be utilized in the hyperbaric center, but the hand sanitizer cannot have a high percentage of alcohol as alcohol-based products cannot go in the chamber.

• Therefore, it is not recommended to have alcohol-based hand sanitizers in the hyperbaric chamber room.







• If the first patient of the day has MRSA, it is adequate to clean the chamber at the end of the day?

TRUE --OR-FALSE



• If the first patient of the day has MRSA, it is adequate to clean the chamber at the end of the day?

FALSE



• Where do you find the list of approved chamber cleaners?



- Where do you find the list of approved chamber cleaners?
 - perrybaromedical.com
 - •sechristusa.com



• How often should the chambers be cleaned?



How often should the chambers be cleaned?

 The chambers should be cleaned with approved cleaners weekly or between patients if the patient has contact precautions due to infection, such as MRSA.



Alcohol-based cleaners are the best for cleaning the mattress.

TRUE -OR-FALSE



Alcohol-based cleaners are the best for cleaning the mattress.

FALSE



- What cloth is best to be used to wipe the acrylic?
 - Paper towels
 - Microfiber towels
 - Cotton fabric




- What cloth is best to be used to wipe the acrylic?
 - Paper towels
 - Microfiber towels
 - Cotton fabric








August 2023 – HBO Treatments






August 2023 — Average Daily Patient Volume 22 Business Days





August 2023 – HBO Show Rate





Think HBO

- Education yourself
- Know all the advanced modalities that are available to your patient
- Utilize HBO Techs to help educate the patients
- Saving Lives, Saving Limbs









MATT SCHWEYER, CPCO, CHT-A, CHWS, UHMSDSA CHIEF QUALITY OFFICER | NATIONAL SAFETY DIRECTOR

C. 214-315-5109 888-960-1343 x1007 mschweyer@serenagroups.com www.serenagroupinc.com 125 Cambridge Park Drive Suite 301 Cambridge MA 02140

